




Gelling Fibre Dressing

- High retention capacity to prevent leakage and maceration
- High tensile strength to enable dressing removal in one piece
- Absorbs and retains exudate, blood and bacteria
- Highly absorbent, also under compression therapy
- Soft and conformable which makes it easy to apply

High retention capacity


- Prevents leakage
- Prevents maceration
- Absorbs and retains exudate, bacteria and blood




PVA Fibre technology

- Excellent tensile strength
- Dressing removal in one piece
- Saves nursing time and increases patient comfort

How to use Exufiber




1. Cleanse the wound with saline solution or water. Dry the surrounding skin thoroughly. Donor sites: Harvest tissue and achieve haemostasis according to normal routine.



2a. Apply a dry Exufiber dressing to the wound.

2b. Loosely pack ribbon or sheet into the wound to allow room for swelling of the dressing.



3. Fixate with an appropriate secondary dressing. Fixate (or cover) with an adhesive dressing or a bandage as appropriate. Compression therapy may be used in conjunction with Exufiber.

How Exufiber works

Exufiber is a sterile nonwoven dressing made from highly absorbent polyvinyl alcohol fibres. In contact with wound exudate, Exufiber transforms into a gel that facilitates moist wound healing and ease of removal during dressing change. Exufiber absorbs and retains wound exudate. Exufiber is available both as sheet and ribbon dressings.

Benefits of Exufiber

- High tensile strength to enable dressing removal in one piece
- High retention capacity to prevent leakage and maceration
- Absorbs and retains exudate, blood and bacteria
- Highly absorbent, also under compression therapy
- Soft and conformable which makes it easy to apply
- Can be left in place for up to 7 days and up to 14 days on donor sites

Precautions

All wounds should be inspected frequently. In case of signs of clinical infection, consult a health care professional for adequate infection treatment.

Exufiber is not intended for dry wounds or full thickness burns.

If the dressing dries out and is difficult to remove, it should be moistened according to local policies (e.g. with sterile saline or sterile water) and allowed to soak until it lifts easily. It may take several minutes for Exufiber to transform into a gel. Remove the dressing by gently cleansing/flushing.

Intended Use

Exufiber wound dressing is intended to be used on a wide range of exuding wounds:

- Leg and foot ulcers
- Pressure injuries
- Partial thickness burns
- Surgical wounds
- Donor sites
- Malignant wounds
- Dermal lesions and other external wounds inflicted by trauma

Configurations available both for flat and cavity wounds



Exufiber assortment (Sterile packed)

Art. no	Size cm	Pcs/shelf cont.	Pcs/transp. cont.
603300	5x5	10	40
603301	10x10	10	40
603302	15x15	10	60
603306	4.5x20	10	50
603303	20x30	5	20
603304	2x50	5	20

Find out more at www.molnlycke.com.au

Mölnlycke Health Care AB, Suite 1.01, 10 Tilley Lane, Frenchs Forest NSW 2086. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231, www.molnlycke.co.nz
The Mölnlycke Health Care, Exufiber® names and respective logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2014 Mölnlycke Health Care AB. All rights reserved. HQWC0011-AU