

# Dressing awkward areas with ALLEVYN<sup>®</sup> Ag

Wounds often present in difficult to dress areas. Providing a patient with a dressing that remains in place despite the awkward position of their wound can be very challenging for a clinician.

Here are some ideas from clinicians around the world on how to help gain control of difficult to dress wounds with ALLEVYN Ag.



## Between toe dressing



Take an ALLEVYN Ag Non Adhesive dressing 5cm x 5cm or 10cm x 10cm.

Fold the dressing in half. Make two incisions into the dressing 1-2cm apart.



Feed the toes into the incisions made.

Ensure there is enough dressing between the toes.



Secure in place with tape.

Ensure tape forms a border to the dressing and does not cover ALLEVYN Ag Non-Adhesive completely.



## Between finger dressing



Take an ALLEVYN Ag Non Adhesive dressing 5cm x 5cm or 10cm x 10cm.

Fold the dressing in half. Make two incisions into the dressing 1-2cm apart.



Feed the fingers into the incisions made. Ensure there is enough dressing between the finger.



Secure in place with tape.

Ensure tape forms a border to the dressing and does not cover ALLEVYN Ag Non-Adhesive completely.

## Rectangular dressing



Take an ALLEVYN Ag Non Adhesive square dressing.

Cut to size.



Secure in place with tape.

Ensure tape forms a border to the dressing and does not cover ALLEVYN Ag Non-Adhesive completely.



## Limb dressing (secured with tubular bandage)



Take an ALLEVYN Ag Non Adhesive Dressing of a suitable size (approx. 2cm wider than the burn).

Apply the dressing to the limb. Secure with a small amount of tape at the top and bottom of the dressings, if necessary.

If extra flexibility is required, make slim incisions into the width of the burn, approximately 4cm deep.


Secure in place with tubular bandage, toe to knee.

## Limb dressing (secured with tape)



Take an ALLEVYN Ag Non Adhesive.

Dressing of a suitable size (approx. 2cm wider than the burn).




Apply the dressing to the limb. If extra flexibility is required, make slim incisions into the dressing the width of the burn, approximately 4cm deep.



Secure in place with tape.

Ensure tape forms a border to the dressing and does not cover ALLEVYN Ag Non-Adhesive completely.



## Ordering Information

### ALLEVYN Ag Non Adhesive

Code	Size	UOM
66800083	5cm x 5cm	10
66800086	10cm x 10cm	10
66800089	15cm x 15cm	10
66800091	20cm x 20cm	10


**Please Note:**  
 Never put tape around the circumference of a limb.

Smith & Nephew Pty Ltd Healthcare Division  
 Australia  
 T 61 3 8540 6777 F 61 3 9544 5086  
 www.smith-nephew.com/australia  
 Customer Service  
 T 13 13 60 F 1800 671 000


\*Trademark of Smith & Nephew SN11467 (06/14)

## Burn requiring dressing larger than 20cm x 20cm

### Option 1



### Option 2



Take two dressings of the same size (i.e. two 15cm x 15cm or two 20cm x 20cm dressings), ensuring that the total size will cover the wound with approx. 2cm margins around the wound.

Either apply dressings separately then apply tape to anchor the two dressings together, or use tape to secure the two dressings together before applying the dressing to the burn.

Ensure tape forms a border to the dressing and does not cover ALLEVYN Ag Non-Adhesive completely.