

CUTICELL® – TRIED AND TESTED QUALITY RANGE

THE CHOICE IS NOW YOURS. THE NEW RANGE OF LOW-ADHERENT
WOUND CONTACT LAYERS

ORDERING INFORMATION

Ref. No.		Size	Items per unit
Cuticell® Sterile and individually sealed			
72539-00		7.5cm x 7.5cm	10/bx
72539-02		7.5cm x 20cm	10/bx
72539-04		20cm x 40cm	25/bx
Cuticell® Classic Sterile and individually sealed			
72538-00		5cm x 5cm	5/bx
72538-01		5cm x 5cm	50/bx
72538-02		10cm x 10cm	10/bx
72538-03		10cm x 10cm	100/bx
72538-06		10cm x 7m	1/tin



Cuticell®

PROTECTING, SOOTHING,
HELPING TO HEAL



BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSN0249 04/2017

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSNmedical

an Essity company

THERAPIES. HAND IN HAND.

BSNmedical

an Essity company

CUTICELL® CLASSIC

THE PARAFFIN-IMPREGNATED DRESSING

Cuticell® Classic is the economical solution when treating exuding wounds. The low-adherent dressing is made of open weave cotton gauze impregnated with soft paraffin – to facilitate dressing change and to support a moist environment for wound healing.

PRODUCT DETAILS

- Impregnated with soft paraffin
- Open weave 100% cotton fabric with interlocking threads

INDICATIONS

- First degree and second degree burns
- Lacerations and abrasions
- Leg ulcers
- Donor and recipient skin graft sites

CUTICELL®

THE OINTMENT DRESSING

Cuticell® is soothing and low-adherent, thus minimising pain at dressing changes. The primary wound contact layer is made of smooth acetate fabric impregnated with non-medicated Eucerit® ointment. The result: an ointment dressing that ensures that the wound is not disturbed and helps maintain a moist environment that supports wound healing.

PRODUCT DETAILS

- Eucerit® based ointment
- Smooth acetate fabric
- Defined small mesh size
- Permeable to air and exudate

INDICATIONS

For superficial exuding wounds such as:

- First and second-degree burns
- Cuts and abrasions
- Lacerations
- Donor and recipient skin graft sites

CUTICELL® OFFERS PRODUCT BENEFITS

HELPS TO MAINTAIN A MOIST ENVIRONMENT TO SUPPORT WOUND HEALING

- Does not provide the wound with excess ointment
- Also suitable for patients with sensitive skin
- Impregnated with non-medicated Eucerit® ointment

PREVENTS TISSUE INGROWTH

- Does not disturb wound healing
- Defined small mesh size of 1 x 1mm

MINIMISES PAIN AT DRESSING CHANGE

- Cuticell® is lint-free and does not fray when cut to size
- Smooth acetate fabric with knit structure


ALLOWS EXUDATE TO DRAIN FREELY INTO AN ABSORBENT SECONDARY DRESSING

- Prevents occlusive effects
- Permeable to air and exudate

CUTICELL® CLASSIC OFFERS RELIABLE TREATMENT QUALITY


SUPPORTS WOUND HEALING

- Low-adherent and soothing
- Facilitates dressing change
- Helps maintain a moist environment to support wound healing
- Impregnated with soft paraffin



INCREASED STABILITY WHEN CUT TO SIZE

- For increased stability when cut to size
- Does not fray
- Interlocking weave gauze fabric



ALLOWS EXUDATE TO DRAIN FREELY INTO AN ABSORBENT SECONDARY DRESSING

- Prevents occlusive effects
- Permeable to air and exudate

