

# DEVILBISS SUCTION UNIT

A Fusion of Quality and Innovation



EBOS Code: 34202145  
Supplier Code: 7305P-I

For more than twenty years, DeVilbiss has been delivering the best portable suction options. The DeVilbiss 7305 Series Portable Suction Unit meets suction and capacity requirements in a small, lightweight design.

The DeVilbiss 7305 series Suction Unit is a compact medical suctioning device which has been designed for reliable, portable operation. Because of the small size, light weight and DC operation, the DeVilbiss 7305 Suction Unit is ideal for providing suction in the clinical setting, or on the go with the DC cord. The 7305P series with internal rechargeable battery is also available.

## PERFORMANCE FEATURES

- High Performance - Vacuum adjustments allow for 50-550 mm Hg and a free flow of 27 litres/minute.
- Adjustable Flow Regulator - Easy-to-turn knob conveniently placed directly below gauge, for easy use and control. Knob has safety lock feature which will not allow the knob to be twisted off.
- The DeVilbiss® Suction Unit 7305P series meets the RTCA/DO160D Aircraft Standard as well as the ISO-10079.1:1991 Suction Standard
- Increased Battery Capabilities – For use up to one hour. This provides the patient with added security when it is needed most.

## DESIGN FEATURES

- Easy-to-Read Gauge - The unit has conveniently placed gauge for easy reading. This is especially important when precise measurements of vacuum are needed for the pediatric patient.
- The 7305 series is equipped with a switch mode power supply allowing operation on any AC voltage (100-240 VAC, 50/60 Hz).
- Two-year warranty on Suction Unit.
- Compact Design - All of the high-performance capabilities are conveniently placed in a small compact package weighing less than 6.3 lbs. The compact design with the highperformance capabilities makes the 7305 series suction unit the smallest and strongest portable suction unit of its kind.
- Low battery and charge indicator lights on portable unit (7305P-D).

## ACCESSORIES:

EBOS CODE:	SUPLIER CODE:	PRODUCT DESCRIPTION
34202145	7305P-I	DeVilbiss Portable H/Performance Suction Unit
34232206	7305D-633	DeVilbiss Collection Disposable Bottle (800ml) and Lid - Kit includes Elbow with Integrated Filter, Connection Tubing and 6' Patient Tubing
34203130	7305D-612	4 in 1 Replacement Tubing pkt 6
34232207	7305D-632S	DeVilbiss Collection Disposable Bottle (800ml) and Lid Kit Includes Elbow with Integrated Filter, Connection Tubing
34203125	6305D-611	DeVilbiss 6' Connecting Patient Tubing
34232208	7305D-608	DeVilbiss VacuAide Bacterial Filter - Pkt 12
34202160	7305P-614	DeVilbiss 12vlt DC High Capacity Battery



FOR MORE INFORMATION PLEASE CONTACT CUSTOMER SERVICES OR YOUR ACCOUNT MANAGER ON **FREEPHONE 0800 18 17 16**

[www.eboshealthcare.co.nz](http://www.eboshealthcare.co.nz)