



# Safety Face Shields: PUTTING YOUR HEALTH PROFESSIONALS FIRST

FACE SHIELDS PROVIDE A BARRIER TO ACUTELY-EXPELLED AEROSOLS OF BODY FLUIDS AND ARE COMMONLY USED AS AN ALTERNATIVE TO GOGGLES AS THEY CONFER PROTECTION TO A LARGER AREA OF THE FACE.

## Disposable Fluid Resistant Full-Face Visor

For the protection of the facial area and associated mucous membranes (eyes, nose, mouth) from splashes, sprays, and spatter of body fluids.

### Features:

- Optically clear
- Wrap around
- Full face visor
- No distortion
- Fluid resistant

### Ideal for:

- Aged Care
- Primary Care



FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER  
ON **FREEPHONE 0800 105 501** OR VISIT [www.eboshealthcare.co.nz](http://www.eboshealthcare.co.nz)

ESSENTIAL TO HEALTH



# Safety Face Shields:

PUTTING YOUR HEALTH PROFESSIONALS FIRST



## Visor Specifications:

- Sized to: 33cm x 22cm
- 2cm elastic Band - domed to each side
- Adhesive Tape on rear to apply Foam
- PET of 0.285 anti fog I Clear
- Supplied with protective film both sides - to be removed

## Foam Specifications:

- Sized to: 23cm x 3.5cm @ 3.3cm depth
- Density: 8-20g / cm<sup>2</sup>
- Colours - White, Blue or Black

**EBOS Code:** 28830001



**EBOS**  
HEALTHCARE

FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER  
ON **FREEPHONE 0800 105 501** OR VISIT [www.eboshealthcare.co.nz](http://www.eboshealthcare.co.nz)

ESSENTIAL TO HEALTH