

Frequently Asked Questions Homecare Suction Unit

Q1: How long will the battery operate the unit?

A1: Approximately 60 minutes.

Q2: What is recommended for a discharged battery?

A2: Fully charge the battery for 10-17 hours (depending on depth of discharge).

Q3: If not in use for an extended period of time are there any recommendations for keeping the battery in good condition?

A3: Yes, it's our recommendation that the battery be recharged a minimum of every 6 months.

Q4: If traveling internationally is there anything needed to use my Homecare Suction Unit?

A4: Yes, the correct power cord must be used, as the unit can operate on any AC voltage.

Q5: How do I get replacement accessories and parts for my Homecare Suction Unit?

A5: Through an Authorized DeVilbiss Provider, which can be found on our web site under Find a Provider, located at the very top of our web pages.

Q6: What is the warranty of my Homecare Suction Unit?

A6: 7305P & 7305D – 2 years from date of purchase; the internal rechargeable battery is warranted for 90 days.

