

# PARI LC® SPRINT BABY Nebuliser

## From the first month of life onwards

### The key features at a glance

- Breath enhanced nebuliser (red nozzle insert) for a very fine droplet spectrum, adapted to the anatomy of babies
- With BABY Mask and BABY Bend for a stress-free inhalation whilst sitting or lying down
- Few components – very quick setup
- Straight forward assembly – prevents incorrect operation
- Practical medication level indicator



**PARI LC® SPRINT BABY Nebuliser** (red nozzle insert) with BABY Bend and BABY mask (4 sizes)

### The aerosol characteristics of the PARI LC® SPRINT BABY Nebuliser in combination with:

#### **PARI BOY® SX**

**Total Output Rate:** 150 mg/min

**MMD:** 2.5 µm  
(mass median diameter)

**Mass percentage below 5 µm:** 82%

Measured with the Malvern MasterSizer X (calculated according to Fraunhofer-model) at 23 °C and 50 % relative humidity. Nebulised medium: 0.9 % NaCl (5 ml). Measured at a baby inspiration flow rate (approx. 6 l/min).

### Versions of the PARI LC® SPRINT BABY Nebuliser

	<p><b>PARI LC® SPRINT BABY Nebuliser:</b> Size 0 (for premature babies) Item No. 023G1400</p>	
	<p><b>PARI LC® SPRINT BABY Nebuliser:</b> Size 1 (age 0 - 1) Item No. 023G1401</p>	
	<p><b>PARI LC® SPRINT BABY Nebuliser:</b> Size 2 (age 1 - 3) Item No. 023G1402</p>	
	<p><b>PARI LC® SPRINT BABY Nebuliser:</b> Size 3 (age 3+) Item No. 023G1403</p>	

The **PARI LC® SPRINT BABY Nebuliser** can be combined with all PARI compressors and the **PARI CENTRAL**.