

TRADITIONAL PEP THERAPY OR COMBINATION PEP/NEBULIZER THERAPY

EBOS Code: 22143078

Designed for use with the PARI LC® Plus or LC® Sprint, combination therapy provides concurrent PEP therapy with clinically proven aerosol delivery. Shorter daily treatment regimens and clinically proven aerosol delivery result in an improved patient quality of life.

FEATURES:

- Eight resistance settings cover a large expiratory flow range that offers pressure of 10-20cm H20
- Dishwasher safe & boilable
- Nose clip, mouth piece and inspriatory valve included
- Without manometer

WHAT IS PEP THERAPY:

- Positive Expiratory Pressure (PEP) is a type of therapy that creates resistance pressure during exhalation
- Used for airway stabilization and mucus clearance
- For patients diagnosed with Cystic Fibrosis or lung diseases with secretory problems.



OPTION ONE: Combined PEP and PARI LC® Reusable Nebulizer Therapy



OPTION TWO: Conventional PEP Therapy

ACCESSORIES:

EBOS CODE:	DESCRIPTION
22143071	PARI LC® Sprint
22143079	PARI Manometer
22149345	PARI Nose Clip

