



The impregnated gauze for superficial, exuding wounds



CUTICERIN° impregnated gauze for superficial, exuding wounds

CUTICERIN is an impregnated, non-aqueous, non-linting dressing ideal as a wound contact layer. It reduces the risk of granulation tissue growing through the dressing, decreasing trauma on removal.

Advantages

Low adherent

Made from a smooth acetate fabric with a defined mesh size and water-repellent ointment coating.

Breathable

CUTICERIN allows air and wound exudate to pass freely through the dressing.

Application

- Place CUTICERIN directly onto the wound and cover with a secondary dressing
- CUTICERIN may be cut to the shape of the wound if required

Presentation

CUTICERIN is sterile and individually sealed.

Indications

For the management of superficial, exuding wounds, such as;

- Minor burns
- · Radiation injuries
- Abrasions
- Split-thickness skin graft donor sites





Ordering information

Code	Product description	Items per unit
66045560	7.5cm x 7.5cm	Box/50
66045561	7.5cm x 20cm	Box/50
66045502	20cm x 40cm	Box/25*

*Not available in New Zealand

Reference: 1. Petres, J. Moist wound dressing in superficial skin defects. Zeuschrift fur Haut Krankhethen 1988; 63 (supp 2): 25-29.

Smith & Nephew Pty Ltd Australia T +61 2 9857 3999 F +61 2 9857 3900 www.smith-nephew.com/australia Customer Service

T 13 13 60 F 1800 671 000

New Zealand
T +64 9 820 2840 F +64 9 820 2841
www.smith-nephew.com/new-zealand
Customer Service
T 0800 657 799 (Surgical) T 0800 807 663 (Wound)

Smith & Nephew Ltd

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.