

# CLEMENT CLARKE INTERNATIONAL

## MINI-WRIGHT AFS LOW RANGE



**MINI-WRIGHT PEAK FLOW METERS ARE FULLY FUNDED ON THE PHARMACEUTICAL SCHEDULE FROM NOVEMBER 1ST 2015**

EBOS Code: 27768506

Pharmacode: 248 9945



HAVING SET THE STANDARD WITH THE MINI-WRIGHT STANDARD PEAK FLOW METER, THE AFS LOW RANGE PROVIDES THE SAME LEVEL OF ACCURACY AND CONSISTENCY TO THOSE PATIENTS WITH LOWER PEAK FLOW MONITORING REQUIREMENTS.

Particularly suitable for very young children, the elderly and patients with severe air flow obstruction the Mini-Wright AFS Low Range measures peak expiratory flow between 30-400 L/min.

Particularly useful for patients who have a typically low peak flow reading. By reducing the range of the scale, this device enables the reading to be more accurate in comparison to a scale on a larger peak flow meter.

### PRODUCT FEATURES AND BENEFITS

- Ideal for children and the elderly
- Incorporates aerodynamic flow systems (AFS)
- Individually calibrated for accuracy
- Easy to read scale
- Simple to use
- Available with ATS scale
- Compact for portability

### PRODUCT SUITABILITY

- Hospital
- Clinic
- Home
- Travel

### TECHNICAL SPECIFICATIONS

Dimensions:	176 x 95 x 40mm
Weight:	96g
Measurement Range:	30-400 L/min (ATS Scale)
Accuracy:	±10% or 10 L/min
Repeatability	< 5 L/min

NB: Now available on the PHARMAC H Schedule.