

MICROTOUCH® Blue Nitrile™

ESSENTIAL SERIES EXAMINATION GLOVES

Nitrile, powder free examination glove with exceptional fit and reliable protection

This non-sterile glove features an advanced stretch formulation providing you with the ultimate comfort and fit. With its high strength and excellent protection features, this glove is ideal for nearly all clinical settings.

KEY FEATURES & BENEFITS

- Advanced formulation for exceptional comfort and fit providing better tactile sensitivity
- An AQL better than industry standards delivers consistent, high quality production
- 200-count packaging reduces waste
- Tested for use with over 55 chemicals and approved for use with chemotherapy drugs*

* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations. Users should test the suitability of this product against specific chemicals and the environment where used. Data for specific chemicals is available upon request.

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

RECOMMENDED FOR

- Environmentally friendly facilities
- Protection from latex Type I allergy in HCW's or patients
- Standard, moderate risk examination procedures
- Tested for use with chemotherapy drugs



MICROTOUCH® Everyday, Everywhere Protection

MICROTOUCH® premium medical examination gloves are scientifically developed to fully equip and comfortably protect medical professionals so that they have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICROTOUCH® Blue Nitrile™

ESSENTIAL SERIES EXAMINATION GLOVES

PRODUCT INFORMATION

Material	Nitrile	Freedom From Holes	1.5 AQL
Not Made From Natural Rubber Latex	Yes	Palm Thickness mm (mils)	0.07 (2.8)
Color	Blue	Allergy Prevention	Latex (Type 1)
Cuff Length	Standard	Tested For Use With Chemotherapy Drugs	Yes in accordance with ASTM D6978
Powder Content	Powder-free		
External Glove Surface	Textured Fingers	Sterile	No

PACKAGING AND STORAGE

Packaging	200 gloves per box; 10 boxes per carton/case
Shelf Life	3 years
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, EN ISO 374-1, EN ISO 374-5, EN 374-2, EN 374-4, EN 455 Part 1-4, ISO 10993-10 :2002, ISO 9001 :2008, ISO 13485 :2003, Personal Protective Equipment: Category III, Medical Device: Class I

